

ISOLATION TECHNOLOGY FOR
**PRIMARY PHARMACEUTICAL
PRODUCTION**



Your primary choice.

ISOLATION TECHNOLOGY FOR **PRIMARY PHARMACEUTICAL PRODUCTION**


Comecer isolators for primary pharmaceutical production have been conceived for pharmaceutical companies who require high containment level during manipulations of active pharmaceutical ingredients (APIs).

**Some types
of our isolators for
APIs handling:**


- Dispensing & Weighing
- Milling & Micronizing
- Reactor Charging
- Powder Processing
- Chemical Synthesis
- R&D
- Kilolabo
- Quality Control

A real tailor made approach


Every single process requires a special attention and a deep understanding of the needs. Therefore the final design is at the end what perfectly fits the customer requirements.




Weighing Isolator




Research & Development Isolator



Reactor
Charging Isolator



Dispensing & Weighing Isolator



Quality Control Isolator

Customized internal solutions

Comecer isolators for primary pharmaceutical production can be equipped with racks, accessories and instruments (scales, QC lab instruments, vacuum ovens etc.) in order to answer any specific requirement.



Continuous Liner



Interface Flange



Reactors

Containment: the heart of APIs production

The toxicity of APIs is measured using the OEB/OEL levels. For each of these levels a different containment system is required based also on the process, volume and dustiness of the APIs. Ultimately the isolation technology is the only solution that can guarantee the protection of the operator even in the case of highly toxic products.

OEB:
Occupational Exposure Band

OEL:
Occupational Exposure Level


PDE:
Permitted Daily Exposure

Isolator Systems	6	<0.1	<1	Extremely Toxic
Isolator Systems or Liner Connecting Systems	5	0.1-1	1-10	Highly Toxic
Isolator Systems or Liner Connecting Systems	4	1-10	10-100	Toxic
Liner Connecting Systems or Protective Liner Connecting System Trickle Protection	3	10-100	100-1000	Less Toxic
Protective Liner Connecting System Trickle Protection	2	100-1000	1000-10000	Almost Nontoxic
Open systems with local extraction	1	1000-5000	>10000	Nontoxic
OEB		OEL [µg/m³]	PDE [µg/day]	


General process cycle

OEB 4/5


1. Material loading/operation




2. Transfer chamber washing (CIP/WIP)



3. Material loading during operation




4. Final washing (CIP/WIP)




OEB 5+


1. Material loading




2. Second liner connection



3. Operation



4. Final washing



W Waste

A API: Active Pharma Ingredient

 Contaminated

 Washed

An example: Modular Dispensing and Weighing Isolator

The Modular Dispensing and Weighing Isolator has been conceived to satisfy the requirements related to the handling process of highly toxic API (Active Pharmaceutical Ingredients).

State of the art features

- Stainless steel AISI 316L enclosure with roughness less than 0,4 micron
- Turbulent airflow or laminar airflow
- Inlet/outlet HEPA 14 safe change filters (bag-in bag-out or push-push)
- Nitrogen blanketing system
- Inflatable gaskets
- Washing in place (WIP) / Cleaning in place (CIP)
- Sanitary drain system
- Software CFR21 part 11/GAMP 5 compliant
- PLC control
- Widescreen HMI with user-friendly interface including USB and ethernet port

Further optional

- Instruments integration
- Automatic glove tester (AGT) onboard
- Rapid transfer port (RTP)
- Continuous liner
- Humidity and temperature sensors
- ATEX compliance
- High containment valves (both active and passive parts)
- Paperless graphic recorder
- Others on request

Technical data

Material	
Shell structure	AISI 316 L
Stainless steel finish	external: Scotch-Brite
	internal: Mirror-Brite
Air classification	
Main chamber	Class C ISO 8 (at rest)
Pre-chamber	Class C ISO 8 (at rest)
Working internal pressure	
Main chamber	-125/-25 Pa (standard: -75 Pa)
Pre-chamber	-125/-25 Pa (standard: -50 Pa)
Filters	
Main chamber	N°1 H14 cartridge filter (inlet) N°2 H14 filters bag-in bag-out (outlet)
Pre-chamber	N°1 H14 cartridge filter (inlet) N°2 H14 filters bag-in bag-out (outlet)
Overall dimensions (w x d x h)	2900 x 850 x 2750 mm
Internal dimensions (w x d x h)	
Main chamber	2000 x 600 x 900 mm
Pre-chamber	800 x 600 x 900 mm
Net weight	1300 kg
Utilities requirements	
Compressed air	6 bar, 25 nL/min
Power supply	230V (1Ph+N+PE) 50/60Hz 32A TN-S
Installed power	2600 W

All data refer to the configuration MODULAR DISPENSING AND WEIGHING ISOLATOR (OEB 4/5)



COMECER S.p.A. - Via Maestri del Lavoro, 90
48014 - Castel Bolognese (RA) - Italy
t: +39 0546 656375 - f: +39 0546 656353
comecer@comecer.com - www.comecer.com



www.comecer.com/isolation