IMILL

CONE MILL



IMILL CONE MILL

THE IMILL SERIES OF CONE MILLS CAN OPERATE EITHER AS FREE-STANDING EQUIPMENT OR FULLY INTEGRATED IN OTHER EQUIPMENT.

The IMILL series of cone mills are suitable for wet or dry products, processed within the pharmaceutical, food, nutraceutical and cosmetic industry.

Typical applications for the IMILL are:

- Size reduction or delumping of raw materials prior to powder mixing or granulation.
- Delumping of wet granules to facilitate the drying step.
- Uniform dry-granule particle size to facilitate compression or capsule filling.



WORKING PRINCIPLES

The solid product is fed into the IMILL by gravity or vacuum; the rotation of the impeller forces the product tangentially through the conical screen without generating any friction.

The product, driven by centrifugal acceleration, passes through the holes of the screen with continuous, fast-flow rate, low fine generation, no overheating and

minimum noise. Afterwards, the product is discharged by gravity and can be collected into bags, drums or conveyed directly to the next process stage.



The required particle reduction size can be achieved by changing the IMILL screen and regulating the speed of the impeller. Conical screens are available in various shapes and mesh sizes to handle size reduction of wet, dry, soft and hard products.

AVAILABLE OPTIONS:

- ATEX VERSION
- TEMPERATURE PROBE
- NOZZLES FOR WASH IN PLACE
- TROLLEY WITH HEIGHT REGULATION

TECHNICAL DATA



| Model | 5 | 8 | 12 |
|--|-----------|-----------|-----------|
| Max capacity (kg/h)* | 400 | 2,000 | 4,000 |
| Motor (kW) | 1.5 | 4.0 | 5.5 |
| Speed range (rpm) | 300-2,600 | 300-2,200 | 300-1,500 |
| A (mm) | 180 | 330 | 500 |
| B (mm) | 132.5 | 290 | 414 |
| C (mm) | 125.8 | 199.7 | 301.5 |
| D (mm) | 222 | 369.5 | 420 |
| *Capacity can vary depending on product characteristics, particle size distribution, screen size and installation context. | | | |



