



De Dietrich
PROCESS SYSTEMS

STATIC NUTSCHE FILTERS

Static liquid-solid separation - that is the purpose of Nutsche filters. Fitted with a filtration media at the bottom and suitable for vacuum or pressure applications, Nutsche filters are a larger scale version of the well known lab Büchner filter

De Dietrich
PROCESS SYSTEMS




Simple liquid-solid separation

Whether manufactured in QVF glass, De Dietrich glass-lined steel, stainless steel or hastelloy, Nutsched filters all follow the same principle. A filter media is placed at the bottom of the vessel, allowing a slurry to be loaded and the solids (usually crystals) to be separated from their mother liquor.



Please contact us to discuss your requirement further and to explore the most suited solution for your application.



Which kind of industry can use this product?

Pharmaceuticals and Fine chemicals

Questions? We are here to help.

If you'd like to talk with a sales representative about purchasing De Dietrich Process Systems's products and services, you can reach us here.



03/08/2024 @ 14:32:44