



De Dietrich
PROCESS SYSTEMS

PHARMA REACTOR

Innovative reactors provide optimal cleanability

De Dietrich
PROCESS SYSTEMS



Advanced solutions for efficient cleaning

CHARACTERISTICS

- A wide range from 63 l up to 630 l
- Optimization of the thermal transfer and the mixing performances (OptiMix®)
- Glass-lined flat cover with excellent surface finish and pad nozzles
- Polished 316L Stainless Steel sheathing on drive and motor
- CIP validation
- Successful test results

ADVANTAGES

- OptiMix design for high agitation performances, more connections and less dead zones
- A range of advanced solutions like:
 - Fused Glass quick release opening (by request)
 - Flat seal outlet valve
 - Pad nozzles
 - Reversed mechanical seal
- High performing spray balls: different materials available (in Alloy or PTFE)

Questions? We are here to help.

If you'd like to talk with a sales representative about purchasing De Dietrich Process Systems's products and services, you can reach us here.



03/08/2024 @ 14:47:43