



Robotic Tray Loader RTL




ROBOTIC TRAY LOADER

Overall Dimensions



RTL



Steriline Robotic Tray Loader RTL is designed to arrange the vials in trays with honeycomb distribution up to 600v/min.

The SCARA robot transport the vials from the pusher into the tray. While the SCARA is transporting the vials, the pusher forms another tray.

The RTL can work both at the end of

the line or with half stoppered vials to be transferred to the freeze drier.

Changeover parts consist only of side guides.

Integrated control system, including PLC and HMI, automatically controls all the process and functional parameters of the machine, according to various glassware sizes and types.

Parameters are stored in different recipes, protected by password, accessible only by supervisor and authorized personnel.

All STERILINE tray loaders are built in compliance with cGMP, GAMP and 21CFR Part11 requirements.





Steriline S.r.l.

Via Tentorio, 30
22100 Como - ITALY
T +39 031.592064
sales@steriline.it
www.steriline.it